

GENERAL LAB EQUIPMENT

SHIMADZU BW-K Series Precision platform/large capacity balances



Features & Benefits

- **Windows® Direct Communication Function**
No software is needed to send data to a PC
- **User-Friendly Low Profile**
Compact UniBloc sensor means lowered platform making it easy to load heavy samples
- **GLP/GMP/ISO9000 Print-out**
Calibration records with date and time can be readily produced
- **Built-in Calibration Weight**
Simple lever operation maintains accuracy
- **Application Functions**
Piece counting, percent display, checkweigh-display, and specific gravity measurement software is included
- **High precision/Large capacity**
Large volumes can be read with a high degree of precision and readability

Polystyrene Weighing Boats



Contamination-free pipetting

- Smooth surface for accurate pour-outs with minimal sample loss
- Rounded corners simplify transfer
- May be used with microwave moisture analyzers
- Resist temperatures to 93°C (200°F)

Description	L x W x D (mm)	Cat. No.
White, small, 500/pkg.	41 x 41 x 8	HW-HS1420A
White, medium, 500/pkg.	89 x 89 x 25	HW-HS1420B
White, large, 500/pkg.	140 x 140 x 22	HW-HS1420C
Blue, small, 500/pkg.	41 x 41 x 8	HW-HS1421A
Blue, medium, 500/pkg.	89 x 89 x 25	HW-HS1421B

Sorvall® RC-6™ Plus Superspeed Centrifuge

Watch your productivity soar!

High Speed, High Capacity Processing

The new Sorvall® RC-6™ Plus superspeed centrifuge is designed for maximum productivity, offering 4 liter-capacity and speeds of up to 22,000 rpm with RCFs to 55,200xg. Fast acceleration and deceleration rates deliver additional time savings.

Versatility

The new Sorvall® RC-6™ Plus superspeed centrifuge is available with an array of



SV-46910

Optimized for Safety

The Sorvall® RC-6™ Plus meets CE and CSA safety standards without the need to anchor the unit to the floor.

For more information, please contact Mandel Technical Support at 1-888-883-3636 or techsupp@mandel.ca

ARE Heating Magnetic Stirrer



ARE Model

Features & Benefits

- Cast aluminum heating plate
Plate has special protective coating to ensure uniform heat distribution over the entire surface and resistance to chemicals
- Efficient design ensures that spills do not reach internal parts

VELP Scientifica's stirring line also includes a Vertex-type digital thermoregulator and laboratory stirrers for use with highly viscous liquids. For more information please contact Mandel Technical Support at 1-888-883-3636 or techsupp@mandel.ca

Multistirrers 6 and 15

Multiposition magnetic stirrers



Features & Benefits

- Magnetic Multistirrer 6 is capable of stirring up to 6 x 85mm diameter vessels
- Multistirrer 15 can stir up to 15 x 64mm diameter vessels
- Multistirrers remain cool even after lengthy continuous operation

DISCONTINUED

risk of tube or bottle leakage.