



TempAssure® PCR Tubes



US-1402-8100



US-1402-2900



US-1402-2800

- Tubes and caps sold separately to accommodate your methods and changing needs
- Choice of dome or flat caps
- Strips with individually attached caps are easy to handle and help avoid cross-contamination
- Mixed colour tubes include blue, orange, red, violet and yellow

DESCRIPTION	CAT. NO.	SPECIAL
Individual 0.2mL tubes with attached dome caps, 1,000/pkg.	US-1402-4300*	\$44.85
Individual 0.2mL tubes with attached flat caps, 1,000/pkg.	US-1402-8100*	45.50
8-strip 0.2mL tubes, separate 8-strip dome caps, 125/pkg.	US-1402-2700*	79.95
8-strip 0.2mL tubes, separate 8-strip flat, optically clear caps, 125/pkg.	US-1402-2500	79.95
8-strip 0.2mL tubes with attached individual dome cap, 120/pkg.	US-1402-2900*	91.65
8-strip 0.2mL tubes with attached 8-strip dome cap, 125/pkg.	US-1402-2800*	91.65

*Also available in mixed colours. To specify mixed colours, replace last digit of Cat. No. with 8.

TempPlate® PCR Plates



TempPlate® I



Seal with TempAssure 12- or 1 2-cap strips, or TempPlate® Film.

DESCRIPTION	CAT. NO.	SPECIAL
TempPlate® I, 96-well, non-skirted, flexible PCR plate with raised rims, natural, 10/pkg.	US-1402-9596*	\$30.55
TempPlate® I, 96-well, non-skirted, flexible PCR plate with raised rims, low profile, natural, 10/pkg.	US-1402-9500	30.55
TempPlate® III, 96-well, semi-skirted, natural, 10/pkg.	US-1402-9700*	30.55
TempPlate® V, 96-well, full skirt, natural, 10/pkg.	US-1402-9800	30.55
TempPlate® II, 96-well, chimney top, natural, 10/pkg.	US-1402-9600*	30.55

*Also available in mixed colours. To specify mixed colours, replace last digit of Cat. No. with 8.

SealRite® Microcentrifuge Tubes

SealRite tubes prevent sample evaporation and accidental tube opening



- 1.5mL tube withstands centrifugation up to 25,000xg
- Tamper evident bags have a sealing strip for safer storage
- Stay closed during heating, freezing, centrifugations and exposure to phenol/chloroform
- Seal securely, yet are easy to open
- Made of polypropylene



DESCRIPTION	QUANTITY	CAT. NO.	SPECIAL
Seal-Rite® 1.5mL tubes			
1.5mL natural tube	500	US-1615-5500*	\$14.95
1.5mL natural tube, sterilized	5 x 100	US-1615-5510	20.15
1.5mL assorted colours in individual bags	5 x 100	US-1615-5599	15.60
Seal-Rite® Locking tubes			
1.5mL locking natural tube	500	US-1415-9100	25.35
0.5mL locking natural tube	1,000	US-1605-4400	50.05


* Also available in individual colours. To specify colour, replace the last digit of the Cat. No. with one of the following: 1=blue, 2=green, 4=red, 5=violet, 6=yellow, 7=amber

Tough-Spots®



- Ultra-thin, polyvinyl labels
- Withstand boiling, autoclaving, freezing, liquid nitrogen
- For temperatures from -196°C to 100°C

Tough-Tags®



- Chemically resistant, pre-cut polyester labels
- Withstand freezing, organic solvents and caustic agents
- For temperatures from -196°C to 80°C

DESCRIPTION	LABEL SIZE	FITS	NO. LABELS	CAT. NO.	PRICE	SPECIAL
Small Tough-Spots	3/8" diameter	0.5mL tubes	1000	US-9185-000X*	\$53.00	\$39.75
Large Tough-Spots	0.5" diameter	1.5/2mL tubes	1000	US-9185-050X*	53.00	39.75
Small Tough-Tags	0.94" L x 0.5" W	0.5mL tubes	1000	US-9188-100X*	54.00	40.50
Large Tough-Tags	1.28" L x 0.5" W	1.5/2mL tubes	1000	US-9187-100X*	53.00	39.75
*To order any of the above, please specify colour by replacing x with one of the following: 0=white, 1=blue, 2=green, 4=red, 6=yellow (3=orange, 5=violet, 7=pink; Tough Spots only)						
Teeny Tough-Tags (white)	0.81" L x 0.28" W	0.2mL tubes	1500	US-9188-8690	62.00	46.50
Sidewall Tough-Tags (white)	1.5" L x 0.25" W	96 well plates	1000	US-9188-7925	62.00	46.50
Laser Flag-Tags® (white)	2.5" L x 0.5" W	0.5, 1.5 and 2mL tubes	1500	US-9135-1000	122.00	91.50

Tough-Tags & Tough-Spots for Laser Printers

- Heat resistant 8½ x 11" sheets stay flat and will not jam in laser printers
- Includes instructions for printing with Microsoft® Word
- Offers the same temperature range and chemical resistance as Tough-Tags and Tough-Spots on rolls

www.mandel.ca
View our complete selection of labels!



DESCRIPTION	LABEL SIZE	FITS	TEMP. RANGE	NO. LABELS	CAT. NO.	PRICE	SPECIAL
TOUGH-TAGS FOR LASER PRINTERS (WHITE POLYESTER)							
Small Tough-Tags	0.94" L x 0.5" W	0.5mL tubes	-40°C to 121°C	2975	US-9185-2240	\$109.00	\$81.75
Large Tough-Tags	1.28" L x 0.5" W	1.5/2mL tubes	-40°C to 121°C	2125	US-9187-2016	109.00	81.75
TOUGH-SPOTS FOR LASER PRINTERS (WHITE OR COLOURS - POLYOLEFIN)							
Small Tough-Spots	3/8" diameter	0.5mL tubes	-196°C to 100°C	3840	US-9185-100X*	109.00	81.75
Large Tough-Spots	0.5" diameter	1.5/2mL tubes	-196°C to 100°C	3840	US-9185-200X*	109.00	81.75

* To order any of the above, please specify colour by replacing x with one of the following: 0=white, 1=blue, 2=green, 3=orange, 4=red, 6=yellow, 8=assorted

Cryo-Babies® & Cryo-Tags®



- Ultra-thin, polyolefin labels, adhere to most plastics, glass & metal
- Withstand boiling, autoclaving, freezing, liquid nitrogen
- For temperatures from -196°C to 150°C
- Withstand soaking at room temperature or 100% relative humidity at 35°C for at least 24 hours

DESCRIPTION	LABEL SIZE	FITS	NO. LABELS	CAT. NO.	PRICE	SPECIAL
Small Laser Cryo-Babies, white, 119/sheet	0.94" L x 0.5" W	0.5mL tubes	2380	US-9187-2380	\$109.00	\$81.75
Cryo-Babies, white	1.28" L x 0.5" W	1.5/2.0mL tubes, std. cryo vials	1000	US-9188-2000	61.00	45.75
Laser Cryo-Babies, white, 85/sheet	1.28" L x 0.5" W	1.5/2.0mL tubes, std. cryo vials	1700	US-9187-1700	109.00	81.75
Cryo-Tags, white	1.5" L x 0.75" W	15 and 50mL tubes	1000	US-9188-1500	61.00	45.75
Laser Cryo-Tags, white, 60/sheet	1.5" L x 0.75" W	15 and 50mL tubes	1200	US-9187-1200	109.00	81.75
Large Laser Cryo-Tags, white, 52/sheet	1.69" L x 0.75" W	Standard cryo vials	1040	US-9187-1100	109.00	81.75

Upright Freezer Drawer Rack

Upright tray rack for standard 2-inch high boxes

- Frame constructed of corrosion-resistant stainless steel
- Trays are lightweight aluminum with easy access handles
- 16-box capacity (4 long and 4 high)
- For standard 5.25 x 5.25 x 2 inch high boxes



DESCRIPTION	CAT. NO.	PRICE	SPECIAL
UTR 442 Upright Freezer Drawer Rack	CBS-UTR442	\$232.00	\$162.95

100-Well Storage Rack Boxes

Durable storage for microtubes and cryogenic vials



- Holds 1.5 to 2.0mL tubes and vials
- Measures 141 x 151 x 57mm
- Pack of 5

DESCRIPTION	CAT. NO.	SPECIAL
Assorted colours (blue, green, purple, yellow and orange)	HW-HS1065	\$30.75

81-Well Storage Boxes

- Polypropylene box fits in standard freezer racks
- Stores 1.5 to 2.0mL microtubes under easy-open, friction-fit lid
- Measures 130 x 130 x 47mm – pack of 5



DESCRIPTION	CAT. NO.	SPECIAL
Assorted colours (blue, green, pink, yellow and orange)	HW-HS23255G	\$27.00

Cardboard Cryogenic Vial Boxes and Partitions




- Boxes can be stored in liquid nitrogen or in mechanical freezers
- Box footprint is 133 x 133mm

DESCRIPTION	HEIGHT	SIZE	CAT. NO.	SPECIAL
Vial Box and Lid	50mm	Pkg. of 12	HW-HS2860C	\$27.00
100-Cell Partition (10 x 10)	-	Pkg. of 12	HW-HS2860E	21.75
81-Cell Partition (9 x 9)	-	Pkg. of 12	HW-HS2860H	21.75

Arctic Squares™ Cryostorage Boxes

- Polycarbonate, holds up to 2mL tubes
- Store from -196 to 121°C for repeated use in mechanical freezers as well as liquid nitrogen!



DESCRIPTION	SIZE	CAT. NO.	SPECIAL
5 x 5 Array Boxes - Red	Pkg. of 8	HW-HS23325A	\$46.50
9 x 9 Array Standard Boxes - Assorted	Pkg. of 4	HW-HS23381B	38.25
10 x 10 Array Boxes - Blue	Pkg. of 4	HW-HS23300D	38.25

True North® Cool Containers

Conveniently nest lid beneath bucket or pan

- Rigid single molded polyurethane
- Available in 4 convenient sizes



DESCRIPTION	COLOUR	CAT. NO.	PRICE
4L Bucket	Black Onyx	HW-HS28714K	\$130.00
1L Pan	Emerald	HW-HS28721G	77.00
4L Pan	Arctic Blue	HW-HS28724B	97.00
9L Pan	Ruby	HW-HS28729R	147.00

Powder-Free Latex Gloves

- Micro-textured satin finish for exceptional wet grip and tactile sensitivity
- High elasticity and tensile strength
- Ambidextrous, 9.5" long with beaded cuff

\$8.35



SIZE	PKG.	QTY.	CAT. NO.
Extra Small	Box	100	US-4910-1100
Small	Box	100	US-4910-2200
Medium	Box	100	US-4910-3300
Large	Box	100	US-4910-4400
Extra Large	Box	100	US-4910-5500

Powder-Free Nitrile Gloves with Aloe

- Aloe vera conditions skin
- 100% latex-free with no natural rubber proteins
- Ambidextrous, 9.5" long with beaded cuff

\$13.95



SIZE	PKG.	QTY.	CAT. NO.
Small	Box	100	US-4918-2200
Medium	Box	100	US-4918-3300
Large	Box	100	US-4918-4400
Extra Large	Box	100	US-4918-5500

Powder-Free Soft Nitrile Gloves

- 100% latex-free with no natural rubber proteins
- Odour-free nitrile
- Ambidextrous, 9.5" long with beaded cuff

\$8.99



SIZE	PKG.	QTY.	CAT. NO.
Extra Small	Box	100	US-4971-1100
Small	Box	100	US-4971-2200
Medium	Box	100	US-4971-3300
Large	Box	100	US-4971-4400
Extra Large	Box	100	US-4971-5500

4-Way Microtube Racks

Link together multiple racks



- Heavy-duty polypropylene racks
- Can hold 4 x 50mL, 12 x 15mL, 32 x 1.5mL or 32 x 0.5mL tubes
- Pack of 5

DESCRIPTION	CAT. NO.	SPECIAL
Assorted colours (blue, green, pink, yellow and orange)	HW-HS29022G	\$30.75

80-Well Microtube Racks



- Easily fill and transfer 80 tubes in compact 5 x 16 arrangement
- Polypropylene rack holds 1.5 to 2.0mL microtubes
- Measures 174 x 95 x 52mm
- Pack of 5

DESCRIPTION	CAT. NO.	SPECIAL
Assorted colours (blue, green, pink, yellow and orange)	HW-HS29025G	\$15.00

96-Well PCR Racks

- Polypropylene rack accommodates individual 0.2mL tubes or strips of 8 or 12 tubes
- Measures 174 x 95 x 52mm
- Pack of 5




DESCRIPTION	CAT. NO.	SPECIAL
Assorted colours (blue, green, pink, yellow and orange)	HW-HS23461A	\$27.00

Arctic Ice PCR Tube Rack

Polypropylene 0.2mL chill rack changes colour from green to yellow as temperature increases. Use to prepare cycle sequencing reactions, restriction digestions, or prevent non-specific primer annealing in 0.2mL tubes or 96-well PCR plates. Measures 5.56" L x 3.9" W x 1.74" H with the locking lid.

Colour changes when the temperature exceeds 7°C!



DESCRIPTION	CAT. NO.	SPECIAL
PCR Tube Rack, yellow/green with natural lid, 2/pkg.	US-2312-2796	\$84.50