

Thermo Scientific Heraeus Centrifuges

Exceptional Performance and Throughput

- Process up to 40 x 50mL conical tubes (40% more) in your swing-out rotor, or up to 14 x 50mL (100% more) in fixed angle configuration, for increased capacity over competitive models in the same bench space.

Unmatched Versatility

- Industry's broadest range of rotors and adapters available, including Thermo Scientific Fiberlite rotors, enabling a minimum of 16 rotor configurations to meet all processing needs.

Outstanding Energy Efficiency

- Up to 40% energy consumption savings on industry standard protocols, such as blood separation or conical tubes processing.

Auto-Lock® III Rotor System

Secure locking system delivers:

- Installation and exchange of rotors with just the push of a button for quick-switch between applications
- Easy chamber access enables cleaning convenience, reducing contamination risk and leading to improved research and a healthier working environment

SMARTSpin® Technology

- Advanced rotor management system optimizes acceleration, braking and residual load imbalance for smooth runs and reproducible separations



Heraeus Multifuge X1R

Tissue Culture Package

- Capacity of 1.6L (4 x 400mL)
- 16 x 50mL tubes up to 5,000rpm (4,696 x g)
- Maximum speed of 15,200rpm (25,830 x g) with fixed angle rotor

PACKAGE SPECIAL
\$7,199
Cat. No. SVC-75518382

Package includes:

- Multifuge X1R Refrigerated Centrifuge
- TX-400 Swinging Bucket Rotor (4 x 400mL)
- 400mL Round Buckets (Set of 4)
- Adapters 4 x 50mL Conical (Set of 4 - Capacity 16)
- Adapters 9 x 15mL Conical (Set of 4 - Capacity 36)



Heraeus Megafuge 40R

High Capacity Tissue Culture Package

- Capacity of 3L (4 x 750mL)
- 40 x 50mL tubes up to 3,500rpm (2,700 x g)
- Maximum speed of 15,200rpm (25,314xg) with fixed angle rotor

PACKAGE SPECIAL
\$8,799
Cat. No. SVC-75197473

Package includes:

- Megafuge 40R Refrigerated Centrifuge
- BIOLiner High Capacity Swinging Bucket Rotor (4 x 750mL)
- BIOLiner Buckets (Set of 4)
- Adapters 10 x 50mL Conical (Set of 4 - Capacity 40)
- Adapters 22 x 15mL Conical (Set of 4 - Capacity 88)



Heraeus Megafuge 40R

Tissue Culture Package

- Capacity of 3L (4 x 750mL)
- 28 x 50mL tubes up to 4,700rpm (4,800 x g)
- Maximum speed of 15,200rpm (25,314xg) with fixed angle rotor

PACKAGE SPECIAL
\$9,499
Cat. No. SVC-75193839

Package includes:

- Megafuge 40R Refrigerated Centrifuge
- TX-750 Swinging Bucket Rotor (4 x 750mL)
- 750mL Round Buckets (Set of 4)
- Adapters 7 x 50mL Conical (Set of 4 - Capacity 28)
- Adapters 14 x 15mL Conical (Set of 4 - Capacity 56)



Heraeus Multifuge X3R

Tissue Culture Package

- Capacity of 3L (4 x 750mL)
- 28 x 50mL tubes up to 4,700rpm (4,800 x g)
- Maximum speed of 15,200rpm (25,314xg) with fixed angle rotor

PACKAGE SPECIAL
\$9,999
Cat. No. SVC-75163839

Package includes:

- Multifuge X3R Refrigerated Centrifuge
- TX-750 Swinging Bucket Rotor (4 x 750mL)
- 750mL Round Buckets (Set of 4)
- Adapters 7 x 50mL Conical (Set of 4 - Capacity 28)
- Adapters 14 x 15mL Conical (Set of 4 - Capacity 56)



Spectrafuge 24D Microcentrifuge

Digital and high speed

- Available in blue, teal, purple, grey or red
- Easy access 24-place rotor
- Powerful, brushless drive
- Optional StripSpin adapter for 0.2mL tubes & strips - FREE with purchase of microcentrifuge
- Speed range: 500rpm to 13,300rpm
- Maximum rcf: 16,300 x g



DESCRIPTION	CAT. NO.	PRICE	SPECIAL
Spectrafuge, 24D grey, complete with 24 x 1.5/2.0mL rotor	NL-C-2400	\$1,956.00	\$1,595.00
Spectrafuge, 24D other colours, complete with 24 x 1.5/2.0mL rotor	NL-C-2400-X*	1,956.00	1,595.00
StripSpin adapter for 0.2mL strip tubes	NL-C-2400-SS	121.00	-

*To specify colour replace "X" with appropriate letter to end of catalogue number: B-blue, T-teal, P-purple, R-red

Hermle Z216-MK Refrigerated Microcentrifuge

- High capacity – 44 place rotor for 1.5 / 2.0mL tubes
- Ten acceleration and deceleration levels
- Set speed (rpm or rcf) in increments of 10
- Maximum rcf: 21,380 x g

With 24-place rotor
(Cat. No. NL-C-0216-MK-24)

\$5,195

With 44-place rotor
(Cat. No. NL-C-0216-MK-44)

\$5,295



Prism High Speed Microcentrifuge



\$1,995

Cat. No. NL-C-2500

- Most powerful on the market
- New low-profile design
- Unique Multi-Flow air cooling technology
- Maximum speed of 15,000rpm /21,200 xg
- Includes 24-place rotor

Hermle Refrigerated Z400K High Capacity Centrifuge

- Maximum speed 13,500rpm/17,310xg
- Automatic soft stop braking
- Microprocessor control with g-force conversion program
- Multiple swing out and angle rotor options

Package includes:

- 4 x 500mL rotor
- Adapters for 15mL tubes (capacity 48)
- Adapters for 50mL tubes (capacity 28)



REFRIGERATED PACKAGE SPECIAL

\$7,199

Cat. No. NL-C-0400-K-P1



MPS 1000 Mini-Plate Spinner

\$595

Cat. No. NL-C-1000

- A compact and economical centrifuge for PCR plates
- Quickly spin down droplets and condensation
- Use before and after thermal cycling to improve PCR yield
- Accepts skirted, non-skirted and all standard PCR plates

Mandel Mini Centrifuge

- Max. speed 6000rpm/2000 x g
- Ideal for quick spins and microfiltration
- Rotors for 6 x 1.5/2.0mL tubes, 2 x 0.2ml strips
- Hinge with stainless steel pin is virtually indestructible



\$299

Cat. No. DEL-1000R



**BUY 1 AND
GET 2ND FOR
\$149.50!**

PROMO CODE: MS-MANDELMINI2

Spectrafuge 6C Compact Research Centrifuge



- High speed, up to 6500rpm/4000xg
- Digital control and accuracy
- Small footprint – 8 x 9"
- Accepts all common vacutainers and 6 x 15mL conical tubes

\$849

Cat. No. NL-C-0060

SealRite® Microcentrifuge Tubes

SealRite tubes prevent sample evaporation and accidental tube opening

- 1.5mL tube withstands centrifugation up to 25,000xg
- Packed in tamper-evident, resealable bags
- Stay closed during heating, freezing, centrifugations and exposure to phenol/chloroform
- Seal securely, yet are easy to open
- Made of polypropylene



**SealRite®
Locking Tube
locking prongs
extend over outer
rim of tubes**

DESCRIPTION	QUANTITY	CAT. NO.	SPECIAL
SEALRITE® 1.5ML TUBES			
1.5mL natural tube	500	US-1615-5500	\$14.95
1.5mL blue tube	500	US-1615-5501*	15.95
1.5mL natural tube, sterilized	5 x 100	US-1615-5510	17.95
1.5mL assorted colours (individual bags)	5 x 100	US-1615-5599	15.95
SEALRITE® LOCKING TUBES			
1.5mL locking natural tube	500	US-1415-9100	19.95
0.5mL locking natural tube	1,000	US-1605-4400	49.95

* Also available in individual colours. To specify colour, replace the last digit of the Cat. No. with one of the following: 1=blue, 2=green, 4=red, 5=violet, 6=yellow, 7=amber