

# Laboratory Labels on Rolls

## Tough-Spots®




- Ultra-thin, polyvinyl labels
- Temperature resistant from -196°C to 80°C
- Withstand boiling, autoclaving, freezing, liquid nitrogen

## Cryo-Babies® & Cryo-Tags®



- Ultra-thin, polyolefin labels, adhere to most plastics, glass & metal
- Temperature resistant from -196°C to 150°C Withstand boiling, autoclaving, freezing, liquid nitrogen
- Withstand soaking at room temperature or 100% relative humidity at 35°C for at least 24 hours

## Tough-Tags®



- Chemically resistant, polyester labels
- Temperature resistant from -40°C to 121°C
- Withstand freezing, organic solvents and caustic agents

DESCRIPTION	LABEL SIZE	FITS	TEMP. RANGE	NO. LABELS	CAT. NO.	SPECIAL
<b>TOUGH-SPOTS POLYVINYL LABELS ON ROLLS (WHITE OR COLOURS)</b>						
Small Tough-Spots	3/8" diameter	0.5mL tubes	-196° to 80°C	1000	HW-9185-000X*	<b>\$31.50</b>
Large Tough-Spots	0.5" diameter	1.5/2mL tubes	-196° to 80°C	1000	HW-9185-050X*	<b>31.50</b>
Tough-Spots for 15mL tubes	0.75" diameter	15mL tubes	-196° to 80°C	500	HW-9185-0750	<b>31.50</b>
Tough-Spots for 50mL tubes	1.0" diameter	50mL tubes	-196° to 80°C	500	HW-9185-0100	<b>31.50</b>
<b>CRYO-BABIES &amp; CRYO-TAGS POLYOLEFIN LABELS ON ROLLS (WHITE)</b>						
Cryo-Babies	1.28" L x 0.5" W	1.5/2.0mL tubes, cryovials	-196°C to 150°C	1000	HW-9188-2000	<b>37.80</b>
Cryo-Tags	1.5" L x 0.75" W	15 and 50mL tubes	-196°C to 150°C	1000	HW-9188-1500	<b>37.80</b>
<b>TOUGH-TAGS POLYESTER LABELS ON ROLLS (WHITE OR COLOURS)</b>						
Small Tough-Tags	0.94" L x 0.5" W	0.5mL tubes	-40°C to 121°C	1000	HW-9188-100X*	<b>31.50</b>
Large Tough-Tags	1.28" L x 0.5" W	1.5/2mL tubes	-40°C to 121°C	1000	HW-9187-100X*	<b>31.50</b>
Teeny Tough-Tags	0.81" L x 0.28" W	0.2mL tubes	-40°C to 121°C	1500	HW-9188-8690	<b>37.80</b>
Sidewall Tough-Tags	1.5" L x 0.25" W	96 well plates	-40°C to 121°C	1000	HW-9188-7925	<b>37.80</b>
Laser Flag-Tags®	2.5" L x 0.5" W	0.5, 1.5 and 2mL tubes	-40°C to 121°C	1500	HW-9135-1000	<b>75.60</b>

\*To order any of the above, please specify colour by replacing X with: 0=white, 1=blue, 2=green, 4=red, 6=yellow, 3=orange, 5=violet, 7=pink, 14=maroon, 16=gold, 17=coffee, 18=silver

# Laboratory Labels on Sheets for Laser Printers



## Laser Tough-Spots® Laser Cryo-Babies® & Cryo-Tags® Laser Tough-Tags®

- Heat resistant 8½ x 11" sheets stay flat and will not jam in laser printers
- Include instructions for printing with Microsoft® Word



DESCRIPTION	LABEL SIZE	FITS	TEMP. RANGE	NO. LABELS	CAT. NO.	SPECIAL
<b>TOUGH-SPOTS POLYOLEFIN LABELS FOR LASER PRINTERS (WHITE OR COLOURS)</b>						
Small Laser Tough-Spots	3/8" diameter	0.5mL tubes	-196°C to 150°C	3840	HW-9185-100X*	<b>\$75.60</b>
Large Laser Tough-Spots	0.5" diameter	1.5/2mL tubes	-196°C to 150°C	3840	HW-9185-200X*	<b>75.60</b>
Laser Tough-Spots for 15mL tubes	0.75" diameter	15mL tubes	-196°C to 150°C	2160	HW-9185-3000	<b>75.60</b>
Laser Tough-Spots for 50mL tubes	1.0" diameter	50mL tubes	-196°C to 150°C	1260	HW-9185-4000	<b>75.60</b>
<b>CRYO-BABIES &amp; CRYO-TAGS POLYOLEFIN LABELS FOR LASER PRINTERS (WHITE)</b>						
Small Laser Cryo-Babies, 119/sheet	0.94" L x 0.5" W	0.5mL tubes	-196°C to 150°C	2380	HW-9187-238X*	<b>75.60</b>
Laser Cryo-Babies, 85/sheet	1.28" L x 0.5" W	1.5/2.0mL tubes, cryovials	-196°C to 150°C	1700	HW-9187-170X*	<b>75.60</b>
Laser Cryo-Tags, 60/sheet	1.5" L x 0.75" W	15 and 50mL tubes	-196°C to 150°C	1200	HW-9187-1200	<b>75.60</b>
Large Laser Cryo-Tags, 52/sheet	1.69" L x 0.75" W	Standard cryovials	-196°C to 150°C	1040	HW-9187-1100	<b>75.60</b>
<b>TOUGH-TAGS POLYESTER LABELS FOR LASER PRINTERS (WHITE)</b>						
Small Laser Tough-Tags	0.94" L x 0.5" W	0.5mL tubes	-40°C to 121°C	2975	HW-9185-224X*	<b>75.60</b>
Large Laser Tough-Tags	1.28" L x 0.5" W	1.5/2mL tubes	-40°C to 121°C	2125	HW-9187-2016	<b>75.60</b>

\* To order any of the above, please specify colour by replacing x with one of the following: 0=white, 1=blue, 2=green, 3=orange, 4=red, 6=yellow, 8=assorted

## Tough-Tags® Station



The portable Tough-Tags® station holds up to six standard size boxes of Tough-Tags®, Cryo-Babies®, and Tough-Spots®. The station also accommodates five Tough-Spots® and Teeny Tough-Spots® boxes.

DESCRIPTION	CAT. NO.	PRICE	SPECIAL
Tough-Tags® Station	HW-HS23460	\$71.00	<b>\$56.80</b>

Simport

# AMPLITUBE™ Thin Wall PCR Tubes



SIM-T325-1

AMPLITUBE™ PCR tubes are made of transparent superior quality grade polypropylene for better viewing of the contents. The ultrathin wall design ensures rapid thermal transfer and a significant reduction in cycle and PCR reaction time.

Certified free of detectable levels of DNase, RNase, pyrogen and DNA-free.

- Strips with individually attached caps are easy to handle and help avoid cross-contamination
- Choice of dome or flat caps



SIM-T320-2

DESCRIPTION	QTY./PACK	CAT. NO.	PRICE	SPECIAL
<b>0.2ML AMPLITUBE 8-TUBE STRIPS WITH ATTACHED CAPS, NATURAL</b>				
8-tube strips with attached dome caps	125	SIM-T320-3N*	\$120.00	<b>\$96.00</b>
8-tube strips with attached flat caps	125	SIM-T320-2N*	120.00	<b>96.00</b>
<b>0.2ML AMPLITUBE INDIVIDUAL PCR TUBES, NATURAL</b>				
Tubes with attached dome caps	1000	SIM-T325-1N*	47.00	<b>37.60</b>
Tubes with attached flat caps	1000	SIM-T325-2N*	47.00	<b>37.60</b>
<b>0.2ML AMPLITUBE REACTION STRIPS, NATURAL</b>				
8-tube strips, tubes only	125	SIM-T320-1N*	74.00	<b>59.20</b>
8-cap strips, dome caps only	125	SIM-T321-1N*	15.00	<b>12.00</b>
8-cap strips, flat caps only.	125	SIM-T321-2N*	15.00	<b>12.00</b>


Also available in colours. Please specify colour by replacing N with: B=blue, G=green, R=red, Y=yellow

# AMPLATE™ Thin Wall PCR Plates

Designed for rapid thermal transfer

## 96-well plates

96-well plates are available with 200µL capacity as standard, in a low profile version with 100µL capacity, reducing the dead space between the sample and the cover, and in a skirted low-profile version with 100µL capacity to provide a bar-coding and labeling area, unavailable in the other plates.



SIM-T323-96LPP

## 384-well plates

384-well plates have 40µL well capacity, designed to be used with reaction volumes from 2 to 30µL. The AMPLATE™ 384-well plate is skirted to allow bar coding on sides for identification, and also to make it compatible with automated fluid handling systems.



SIM-T323-384SKN

DESCRIPTION	WELL CAPACITY	QTY./PACK	CAT. NO.	PRICE	SPECIAL
<b>AMPLATE THIN WALL PCR PLATES, NATURAL</b>					
96-well thin wall PCR plates	200µL	100	SIM-T323-96N*	\$265.00	<b>\$212.00</b>
96-well thin wall PCR plates, low profile	100µL	100	SIM-T323-96LPN*	265.00	<b>212.00</b>
96-well thin wall PCR plates, skirted	100µL	100	SIM-T323-96SKN*	306.00	<b>244.80</b>
384-well thin wall PCR plates, skirted	40µL	100	SIM-T323-384SKN*	530.00	<b>424.00</b>

Also available in colours. Please specify colour by replacing N with: B=blue, G=green, R=pink, Y=yellow

Certified free of detectable levels of DNase, RNase, pyrogen and DNA-free.

# Parafilm® Laboratory Sealing Film

Colourless, moisture-proof film moulds to surfaces

Stretch film tight to seal beakers, flasks, tubes, and Petri dishes.

Parafilm is a registered trademark of American Can Company.

DESCRIPTION	CAT. NO.	PRICE	SPECIAL
Parafilm® 50mm x 75m	HW-HS234526A	\$37.00	<b>\$22.20</b>
Parafilm® 100mm x 38m	HW-HS234526B	38.00	<b>22.80</b>
Parafilm® 100mm x 75m	HW-HS234526C	66.00	<b>39.60</b>



# Stainless Steel Freezer Racks



US-2632-4400

US-2632-3300



US-2602-4500

US-2602-4500

US-2602-4400


US-2602-3400

US-2602-4400

### Racks for Chest Freezers

CAT. NO.	DESCRIPTION	PRICE	SPECIAL
<b>Standard 2" Box* Racks</b>			
US-2602-3300	3 x 3, 16 <sup>5</sup> / <sub>16</sub> " x 5 <sup>1</sup> / <sub>2</sub> " x 6 <sup>5</sup> / <sub>16</sub> " H	\$139.00	<b>\$83.40</b>
US-2602-3400	3 x 4, 16 <sup>5</sup> / <sub>16</sub> " x 5 <sup>1</sup> / <sub>2</sub> " x 8 <sup>7</sup> / <sub>16</sub> " H	168.00	<b>100.80</b>
US-2602-3500	3 x 5, 16 <sup>5</sup> / <sub>16</sub> " x 5 <sup>1</sup> / <sub>2</sub> " x 11" H	177.00	<b>106.20</b>
US-2602-4400	4 x 4, 21 <sup>5</sup> / <sub>16</sub> " x 5 <sup>1</sup> / <sub>2</sub> " x 8 <sup>7</sup> / <sub>16</sub> " H	224.00	<b>134.40</b>
US-2602-4500	4 x 5, 21 <sup>5</sup> / <sub>16</sub> " x 5 <sup>1</sup> / <sub>2</sub> " x 11" H	291.00	<b>174.60</b>
<b>2" Box* Racks with Drawers – Load by sliding boxes through the sides of the frame</b>			
US-2632-3300	3 x 3, 16 <sup>1</sup> / <sub>2</sub> " x 5 <sup>1</sup> / <sub>2</sub> " x 7 <sup>1</sup> / <sub>8</sub> " H	183.00	<b>109.80</b>
US-2632-3400	3 x 4, 16 <sup>1</sup> / <sub>2</sub> " x 5 <sup>1</sup> / <sub>2</sub> " x 9 <sup>7</sup> / <sub>16</sub> " H	245.00	<b>147.00</b>
US-2632-4300	4 x 3, 22" x 5 <sup>1</sup> / <sub>2</sub> " x 7 <sup>1</sup> / <sub>8</sub> " H	245.00	<b>147.00</b>
US-2632-4400	4 x 4, 22" x 5 <sup>1</sup> / <sub>2</sub> " x 9 <sup>7</sup> / <sub>16</sub> " H	329.00	<b>197.40</b>
<b>Top-Load 2" Box* Racks with Drawers – Load boxes from the top instead of the side</b>			
US-2687-4400	4 x 4, 22.34" x 5.47" x 9.29" H	239.00	<b>143.40</b>

\*Cardboard or polypropylene boxes measuring 5<sup>1</sup>/<sub>4</sub>" L x 5<sup>1</sup>/<sub>4</sub>" W x 2" H boxes.



US-2602-1100

### Racks for Upright Freezers

STANDARD BOX* RACKS WITH LOCKING RODS	DESCRIPTION	PRICE	SPECIAL
US-2602-1100	1 x 11, 5 <sup>1</sup> / <sub>2</sub> " x 5 <sup>5</sup> / <sub>8</sub> " x 24 <sup>1</sup> / <sub>4</sub> " H	\$154.00	<b>\$92.40</b>
US-2602-1200	1 x 12, 5 <sup>1</sup> / <sub>2</sub> " x 5 <sup>5</sup> / <sub>8</sub> " x 26 <sup>3</sup> / <sub>8</sub> " H	168.00	<b>100.80</b>

\*Cardboard or polypropylene boxes measuring 5<sup>1</sup>/<sub>4</sub>" L x 5<sup>1</sup>/<sub>4</sub>" W x 2" H boxes.

## Cardboard Cryogenic Vial Boxes

- Boxes can be stored in liquid nitrogen or in mechanical freezers
- 133 (W) x 133 (D) x 50 (H)mm



DESCRIPTION	CAT. NO.	SPECIAL
Vial Box and Lid, pk./12	HW-HS2860C	<b>\$23.40</b>
100-Cell Partition (10 x 10), pk./12	HW-HS2860E	<b>18.85</b>
81-Cell Partition (9 x 9), pk./12	HW-HS2860H	<b>18.85</b>

## 15 & 50 mL Storage Box

- Rugged polypropylene construction
- Both 15 and 50mL tube grids included with each box



DESCRIPTION	CAT. NO.	SPECIAL
15 & 50mL Storage Box, pk./5	HW-HS23271A	<b>\$39.00</b>

## Arctic Squares™ Cryostorage Boxes

- Polycarbonate, holds up to 2mL tubes
- Store from -196 to 121°C for repeated use in mechanical freezers as well as liquid nitrogen!



DESCRIPTION	CAT. NO.	SPECIAL
5 x 5 Array Boxes – Red, pk./8	HW-HS23325A	<b>\$40.30</b>
9 x 9 Array Boxes – Assorted, pk./4	HW-HS23381B	<b>33.15</b>
10 x 10 Array Boxes – Blue, pk./4	HW-HS23300D	<b>33.15</b>

## True North® Cool Containers

- Rigid single moulded polyurethane
- Available in 4 convenient sizes



DESCRIPTION	COLOUR	CAT. NO.	SPECIAL
4L Bucket	Black Onyx	HW-HS28714K	<b>\$42.35</b>
1L Pan	Emerald	HW-HS28721G	<b>53.35</b>
4L Pan	Arctic Blue	HW-HS28724B	<b>80.85</b>
9L Pan	Ruby	HW-HS28729R	<b>109.45</b>

# Cryovial® Tubes


For storing biological material as low as -196°C

- Polypropylene, sterilized by gamma radiation
- Packaged in tamperproof, resealable, safety-lock bags of 100
- Should be used only in gas phase of liquid nitrogen



VOLUME (mL)	SIZE (mm)	SELF-STANDING	ROUND BOTTOM	QTY./CS.	CAT. NO.	PRICE	SPECIAL
<b>EXTERNAL THREADS</b>							
1.2	12.5 x 42	•		1,000	SIM-T310-1A	\$346.00	<b>\$276.80</b>
2	12.5 x 47		•	1,000	SIM-T310-2	346.00	<b>276.80</b>
2	12.5 x 49	•		1,000	SIM-T310-2A	346.00	<b>276.80</b>
3	12.5 x 71	•		1,000	SIM-T310-3A	350.00	<b>280.00</b>
4	12.5 x 77	•		1,000	SIM-T310-4A	387.00	<b>309.60</b>
5	12.5 x 91	•		1,000	SIM-T310-5A	423.00	<b>338.40</b>
10	17 x 84	•		500	SIM-T310-10A	290.00	<b>232.00</b>
<b>INTERNAL THREADS</b>							
1.2	12.5 x 41	•		1,000	SIM-T311-1	\$333.00	<b>\$266.40</b>
2	12.5 x 49	•		1,000	SIM-T311-2	333.00	<b>266.40</b>
2	12.5 x 48		•	1,000	SIM-T311-3	333.00	<b>266.40</b>
4	12.5 x 70		•	1,000	SIM-T311-4	382.00	<b>305.60</b>
4	12.5 x 72	•		1,000	SIM-T311-4A	382.00	<b>305.60</b>
5	12.5 x 90		•	1,000	SIM-T311-5	444.00	<b>355.20</b>

1. A Capinsert™ is available in 11 different colours. Perfect for colour coding.
2. Vertical ribs make cap easy to remove.
3. Super fast 1¼ turn thread design.
4. Made of same material, therefore same coefficient of expansion ensures secure seal at all temperatures.
5. Thick wall makes vial almost unbreakable.
6. Large white marking area.
7. Excellent clarity makes sample easy to see.
8. Round bottom. Very easy to empty contents completely.
9. Many sizes available as self-standing with universal locking base.



# Micrewtubes®

Model T335 with o-ring seal and flat screw cap for use down to -196°C

O-ring seal screw caps are unattached but incorporate a flat top which provides a writing surface. These flat surfaced caps can be used with automatic capping machines in packaging industries. Closures are supplied in natural color. Perfect for cryogenic work. Non skirted tubes can be centrifuged up to 20,000g. Will withstand temperatures from -196 °C to +121 °C. With o-ring seal and flat top.



SELF-STANDING 0.5ML	SELF-STANDING 1.5ML	WITHOUT SKIRT 1.5ML	SELF-STANDING 1.8ML	WITHOUT SKIRT 1.8ML	QTY./CS.	PRICE	SPECIAL
Caps not assembled. Non-sterile, non-graduated.							
SIM-T335-2	SIM-T335-4	SIM-T335-5	SIM-T335-6	SIM-T335-7	1,000	\$159.00	<b>\$127.20</b>

## 100-Well Storage Rack Boxes

- Holds 1.5 to 2.0mL tubes and vials
- Pack of 5
- NEW! Black colour!



DESCRIPTION	CAT. NO.	SPECIAL
Assorted colours (blue, green, purple, yellow and orange), pk./5	HW-HS1065	<b>\$27.30</b>
Black, pk./5	HW-HS1063	<b>26.95</b>

## 81-Well Storage Boxes



- Polypropylene box fits in standard freezer racks
- Stores 1.5 to 2.0mL microtubes under easy-open, friction-fit lid
- Measures 130 x 130 x 47mm



DESCRIPTION	CAT. NO.	SPECIAL
Assorted colours (blue, green, pink, yellow and orange), pk./5	HW-HS23255G	<b>\$23.40</b>
Black, pk./5	HW-HS23255H	<b>23.65</b>

## 4-Way Microtube Racks



- Can hold 4 x 50mL, 12 x 15mL, 32 x 1.5mL or 32 x 0.5mL tubes
- Pack of 5

DESCRIPTION	CAT. NO.	SPECIAL
Assorted colours (blue, green, pink, yellow and orange)	HW-HS29022G	<b>\$26.65</b>

## 80-Well Microtube Racks


- Easily fill and transfer 80 tubes in compact 5 x 16 arrangement
- Polypropylene rack holds 1.5 to 2.0mL microtubes
- Pack of 5



DESCRIPTION	CAT. NO.	SPECIAL
Assorted colours (blue, green, pink, yellow and orange)	HW-HS29025G	<b>\$13.00</b>

## Arctic Ice PCR Tube Rack

Colour changes when the temperature exceeds 7°C!



DESCRIPTION	CAT. NO.	SPECIAL
PCR Tube Rack, yellow/green with natural lid, 2/pkg.	US-2312-2796	<b>\$84.50</b>


## 96-Well PCR Racks

- Polypropylene rack accommodates individual 0.2mL tubes or strips of 8 or 12 tubes
- Pack of 5



DESCRIPTION	CAT. NO.	SPECIAL
Assorted colours (blue, green, pink, yellow and orange)	HW-HS23461A	<b>\$23.40</b>

## Microscope Slide Boxes<sup>©</sup>



- Cork interior
- Colour-coded for easy identification

DESCRIPTION	CAT. NO.	SPECIAL
25-Place Slide Box, Blue	HW-HS15989A	<b>\$4.50</b>
50-Place Slide Box, Blue	HW-HS15996A	<b>5.00</b>
100-Place Slide Box, Blue	HW-HS15994A	<b>5.50</b>

## Disposable Staining Slide Trays



- Tray holds 8 slides
- Pack of 4 trays and 1 lid

DESCRIPTION	CAT. NO.	SPECIAL
Disposable Slide Trays	HW-HS15951A	<b>\$32.50</b>

## OneRack<sup>®</sup> Test Tube Racks

- No assembly required and no detachable pieces
- Delrin<sup>®</sup> fiberglass-reinforced polyoxymethylene (POM) for exceptional chemical and heat resistance
- Delrin<sup>®</sup> racks are autoclavable
- Case of 8



ONERACK<sup>®</sup> TEST TUBE RACKS – BLUE

DESCRIPTION	CAT. NO.	SPECIAL
13mm (72 places), 104 x 202 x 59mm	HW-HS27511B	<b>\$112.50</b>
16mm (72 places), 127 x 250 x 70mm	HW-HS27512B	<b>112.50</b>
20mm (40 places), 100 x 252 x 83mm	HW-HS27513B	<b>112.50</b>
25mm (40 places), 120 x 300 x 92mm	HW-HS27514B	<b>122.25</b>
30mm (24 places), 110 x 282 x 85mm	HW-HS27515B	<b>122.25</b>



EXCLUSIVE TO MANDEL

EXCLUSIVE TO MANDEL

Shentex facilities fully comply with cGMP and QSR regulations. They are ISO 9001, ISO 13485, CE certified, FDA approved, and meet AQL 1.5 standards.

Exceptional barrier protection, while maintaining flexibility and tactility, these powder-free gloves are highly tear and chemical resistant. Nitrile gloves are available in a variety of coatings and surface finishes.

## pH Natural® 5.5 Powder-Free Soft Nitrile Exam Gloves

**\$7.50**

- pH 5.5 interior lining moisturizes dry, chapped hands
- Odour-free nitrile
- Dark blue, ambidextrous, 9.5" long with beaded cuff

SIZE	PKG.	QTY.	CAT. NO.
Extra Small	Box	100	STX-N125331/1
Small	Box	100	STX-N125332/1
Medium	Box	100	STX-N125333/1
Large	Box	100	STX-N125334/1
Extra Large	Box	100	STX-N125335/1



## AloeForm® Powder-Free Soft Nitrile Exam Gloves with Aloe Vera

**\$8.95**

- Ultimate tactile sensitivity experience
- Moisturizing ACTIValoe® interior lining
- Green, ambidextrous, 9.5" long with beaded cuff

SIZE	PKG.	QTY.	CAT. NO.
Extra Small	Box	100	STX-N128401/1
Small	Box	100	STX-N127402/1
Medium	Box	100	STX-N127403/1
Large	Box	100	STX-N127404/1
Extra Large	Box	100	STX-N127405/1



**PURCHASE 10 BOXES OF GLOVES, RECEIVE**

**\$5 Starbucks Card**

PROMO CODE: MS-STAR5

## Pearlescent® Powder-Free Glossy Blue Nitrile Exam Gloves

**\$10.95**

- Glossy, Pearlescent® outer finish
- Super soft interior lining for added comfort
- Pearlescent® Blue, ambidextrous, 9.5" long with beaded cuff

SIZE	PKG.	QTY.	CAT. NO.
Extra Small	Box	100	STX-N104341/1
Small	Box	100	STX-N104342/1
Medium	Box	100	STX-N104343/1
Large	Box	100	STX-N104344/1
Extra Large	Box	100	STX-N104345/1



**PURCHASE 10 BOXES OF GLOVES, RECEIVE**

**\$10 Starbucks Card**

PROMO CODE: MS-STAR10

## PowerForm® 12" Extra-Long Powder-Free Exam Gloves

**\$8.95**

- Extended cuff - ideal for high-risk work environments where greater precaution is necessary
- Superior comfort and reliability
- Teal, ambidextrous, 12" long with beaded cuff

SIZE	PKG.	QTY.	CAT. NO.
Extra Small	Box	50	STX-N127361/1
Small	Box	50	STX-N127362/1
Medium	Box	50	STX-N127363/1
Large	Box	50	STX-N127364/1
Extra Large	Box	50	STX-N127365/1



**PURCHASE 10 BOXES OF GLOVES, RECEIVE**

**\$5 Starbucks Card**

PROMO CODE: MS-STAR5