

Powder-Free Soft Nitrile Gloves

- 100% latex-free with no natural rubber proteins
- Odour-free nitrile
- Ambidextrous, 9.5" long with beaded cuff

SIZE	PKG.	QTY.	CAT. NO.
Extra Small	Box	100	US-4971-1100
Small	Box	100	US-4971-2200
Medium	Box	100	US-4971-3300
Large	Box	100	US-4971-4400
Extra Large	Box	100	US-4971-5500

\$8.99



Powder-Free Nitrile Gloves with Aloe

- Aloe vera conditions skin
- 100% latex-free with no natural rubber proteins
- Ambidextrous, 9.5" long with beaded cuff

SIZE	PKG.	QTY.	CAT. NO.
Small	Box	100	US-4918-2200
Medium	Box	100	US-4918-3300
Large	Box	100	US-4918-4400
Extra Large	Box	100	US-4918-5500

\$13.95



Powder-Free Latex Gloves

- Micro-textured satin finish for exceptional wet grip and tactile sensitivity
- High elasticity and tensile strength
- Ambidextrous, 9.5" long with beaded cuff

SIZE	PKG.	QTY.	CAT. NO.
Extra Small	Box	100	US-4910-1100
Small	Box	100	US-4910-2200
Medium	Box	100	US-4910-3300
Large	Box	100	US-4910-4400
Extra Large	Box	100	US-4910-5500

\$8.35



SealRite® Microcentrifuge Tubes

SealRite tubes prevent sample evaporation and accidental tube opening



- 1.5mL tube withstands centrifugation up to 25,000xg
- Packed in tamper-evident, resealable bags
- Stay closed during heating, freezing, centrifugations and exposure to phenol/chloroform
- Seal securely, yet are easy to open
- Made of polypropylene

\$12.95

DESCRIPTION	QUANTITY	CAT. NO.	SPECIAL
Seal-Rite® 1.5mL natural tubes	500	US-1615-5500*	\$12.95
Seal-Rite® 1.5mL natural tubes, sterilized	5 x 100	US-1615-5510	16.95

* Also available in individual colours. To specify colour, replace the last digit of the Cat. No. with one of the following: 1=blue, 2=green, 4=red, 5=violet, 6=yellow, 7=amber


Tough-Spots®



- Ultra-thin, polyvinyl labels
- Withstand boiling, autoclaving, freezing, liquid nitrogen
- For temperatures from -196°C to 100°C



Tough-Tags®



- Chemically resistant, pre-cut polyester labels
- Withstand freezing, organic solvents and caustic agents
- For temperatures from -40°C to 121°C

DESCRIPTION	LABEL SIZE	FITS	NO. LABELS	CAT. NO.	PRICE	SPECIAL
Small Tough-Spots	3/8" diameter	0.5mL tubes	1000	US-9185-000X*	\$53.00	\$37.10
Large Tough-Spots	0.5" diameter	1.5/2mL tubes	1000	US-9185-050X*	53.00	37.10
Small Tough-Tags	0.94" L x 0.5" W	0.5mL tubes	1000	US-9188-100X*	54.00	37.80
Large Tough-Tags	1.28" L x 0.5" W	1.5/2mL tubes	1000	US-9187-100X*	53.00	37.10
*To order any of the above, please specify colour by replacing x with one of the following: 0=white, 1=blue, 2=green, 4=red, 6=yellow (3=orange, 5=violet, 7=pink; Tough Spots only)						
Teeny Tough-Tags (white)	0.81" L x 0.28" W	0.2mL tubes	1500	US-9188-8690	62.00	43.40
Sidewall Tough-Tags (white)	1.5" L x 0.25" W	96 well plates	1000	US-9188-7925	62.00	43.40
Laser Flag-Tags® (white)	2.5" L x 0.5" W	0.5, 1.5 and 2mL tubes	1500	US-9135-1000	122.00	85.40

Tough-Tags & Tough-Spots for Laser Printers

- Heat resistant 8½ x 11" sheets stay flat and will not jam in laser printers
- Includes instructions for printing with Microsoft® Word
- Offers the same temperature range and chemical resistance as Tough-Tags and Tough-Spots on rolls

www.mandel.ca
View our complete selection of labels!



DESCRIPTION	LABEL SIZE	FITS	TEMP. RANGE	NO. LABELS	CAT. NO.	PRICE	SPECIAL
TOUGH-TAGS FOR LASER PRINTERS (WHITE POLYESTER)							
Small Tough-Tags	0.94" L x 0.5" W	0.5mL tubes	-40°C to 121°C	2975	US-9185-2240	\$109.00	\$81.75
Large Tough-Tags	1.28" L x 0.5" W	1.5/2mL tubes	-40°C to 121°C	2125	US-9187-2016	109.00	81.75
TOUGH-SPOTS FOR LASER PRINTERS (WHITE OR COLOURS - POLYOLEFIN)							
Small Tough-Spots	3/8" diameter	0.5mL tubes	-196°C to 100°C	3840	US-9185-100X*	109.00	81.75
Large Tough-Spots	0.5" diameter	1.5/2mL tubes	-196°C to 100°C	3840	US-9185-200X*	109.00	81.75

* To order any of the above, please specify colour by replacing x with one of the following: 0=white, 1=blue, 2=green, 3=orange, 4=red, 6=yellow, 8=assorted

Cryo-Babies® & Cryo-Tags®



- Ultra-thin, polyolefin labels, adhere to most plastics, glass & metal
- Withstand boiling, autoclaving, freezing, liquid nitrogen
- For temperatures from -196°C to 150°C
- Withstand soaking at room temperature or 100% relative humidity at 35°C for at least 24 hours

DESCRIPTION	LABEL SIZE	FITS	NO. LABELS	CAT. NO.	PRICE	SPECIAL
Small Laser Cryo-Babies, white, 119/sheet	0.94" L x 0.5" W	0.5mL tubes	2380	US-9187-2380	\$109.00	\$81.75
Cryo-Babies, white	1.28" L x 0.5" W	1.5/2.0mL tubes, std. cryo vials	1000	US-9188-2000	61.00	45.75
Laser Cryo-Babies, white, 85/sheet	1.28" L x 0.5" W	1.5/2.0mL tubes, std. cryo vials	1700	US-9187-1700	109.00	81.75
Cryo-Tags, white	1.5" L x 0.75" W	15 and 50mL tubes	1000	US-9188-1500	61.00	45.75
Laser Cryo-Tags, white, 60/sheet	1.5" L x 0.75" W	15 and 50mL tubes	1200	US-9187-1200	109.00	81.75
Large Laser Cryo-Tags, white, 52/sheet	1.69" L x 0.75" W	Standard cryo vials	1040	US-9187-1100	109.00	81.75

Trio™ Pipette Stand



- Thick body material and a sturdy, reinforced intermediate plate
- Suitable for all current models of Gilson pipettes (except Pipetman Concept)

Acrylic Pipette Holder



- Upright storage for prolonged life of pipettes
- Saves space – dimensions: 10.6 x 5.9 x 5.9"
- Made of durable clear acrylic
- Accepts numerous pipette makes

Single™ Pipette Holder



- Adheres strongly to any clean and dry surface
- Attach to a rough surface or wall using screws
- Accommodates Gilson Microman®, Distriman® and Pipetman® pipettes (except Pipetman Concept®)
- No adapter clip required

DESCRIPTION	CAT. NO.	PRICE	BUY 3 OR MORE
Single pipette holder, 3/pkg.	MS-61406	\$45.00	\$39.00
Trio, portable rack for 3 pipettes	GF-F161405	92.00	86.85
Acrylic pipette holder, each	MS-23810	49.00	45.00

100-Well Storage Rack Boxes

- Holds 1.5 to 2.0mL tubes and vials
- Pack of 5
- NEW! Black colour!



DESCRIPTION	CAT. NO.	SPECIAL
Assorted colours (blue, green, purple, yellow and orange), pk./5	HW-HS1065	\$27.30
Black, pk./5	HW-HS1063	26.95

81-Well Storage Boxes

- Polypropylene box fits in standard freezer racks
- Stores 1.5 to 2.0mL microtubes under easy-open, friction-fit lid
- Measures 130 x 130 x 47mm



DESCRIPTION	CAT. NO.	SPECIAL
Assorted colours (blue, green, pink, yellow and orange), pk./5	HW-HS23255G	\$23.40
Black, pk./5	HW-HS23255H	23.65

Cardboard Cryogenic Vial Boxes

- Boxes can be stored in liquid nitrogen or in mechanical freezers
- 133 (W) x 133 (D) x 50 (H)mm



DESCRIPTION	CAT. NO.	SPECIAL
Vial Box and Lid, pk./12	HW-HS2860C	\$23.40
100-Cell Partition (10 x 10), pk./12	HW-HS2860E	18.85
81-Cell Partition (9 x 9), pk./12	HW-HS2860H	18.85

15 & 50 mL Storage Box


- Rugged polypropylene construction
- Both 15 and 50mL tube grids included with each box



DESCRIPTION	CAT. NO.	SPECIAL
15 & 50mL Storage Box, pk./5	HW-HS23271A	\$39.00

Arctic Squares™ Cryostorage Boxes

- Polycarbonate, holds up to 2mL tubes
- Store from -196 to 121°C for repeated use in mechanical freezers as well as liquid nitrogen!



DESCRIPTION	CAT. NO.	SPECIAL
5 x 5 Array Boxes – Red, pk./8	HW-HS23325A	\$40.30
9 x 9 Array Boxes – Assorted, pk./4	HW-HS23381B	33.15
10 x 10 Array Boxes – Blue, pk./4	HW-HS23300D	33.15

True North® Cool Containers

- Rigid single moulded polyurethane
- Available in 4 convenient sizes



DESCRIPTION	COLOUR	CAT. NO.	SPECIAL
4L Bucket	Black Onyx	HW-HS28714K	\$42.35
1L Pan	Emerald	HW-HS28721G	53.35
4L Pan	Arctic Blue	HW-HS28724B	80.85
9L Pan	Ruby	HW-HS28729R	109.45

Racks with Drawers

- Reach boxes at back of freezer without removing entire rack
- Holds one 5.25 x 5.25 x 2" box per section
- Corrosion resistant stainless steel
- Boxes sold separately



4 x 4

3 x 3

SIZE	DIMENSIONS (L x W x H)	CAT. NO.	SPECIAL
3 x 3	16.5 x 5.5 x 7.125"	US-2632-3300	\$109.80
3 x 4	16.5 x 5.5 x 9.4375"	US-2632-3400	147.00
4 x 3	22 x 5.5 x 7.125"	US-2632-4300	147.00
4 x 4	22 x 5.5 x 9.4375"	US-2632-4400	197.40

4-Way Microtube Racks

- Can hold 4 x 50mL, 12 x 15mL, 32 x 1.5mL or 32 x 0.5mL tubes
- Pack of 5



DESCRIPTION	CAT. NO.	SPECIAL
Assorted colours (blue, green, pink, yellow and orange)	HW-HS29022G	\$26.65

80-Well Microtube Racks

- Easily fill and transfer 80 tubes in compact 5 x 16 arrangement
- Polypropylene rack holds 1.5 to 2.0mL microtubes
- Pack of 5



DESCRIPTION	CAT. NO.	SPECIAL
Assorted colours (blue, green, pink, yellow and orange)	HW-HS29025G	\$13.00

96-Well PCR Racks

- Polypropylene rack accommodates individual 0.2mL tubes or strips of 8 or 12 tubes
- Pack of 5




DESCRIPTION	CAT. NO.	SPECIAL
Assorted colours (blue, green, pink, yellow and orange)	HW-HS23461A	\$23.40

Arctic Ice PCR Tube Rack

DESCRIPTION	CAT. NO.	SPECIAL
PCR Tube Rack, yellow/green with natural lid, 2/pkg.	US-2312-2796	\$84.50

Colour changes when the temperature exceeds 7°C!



Disposable Staining Slide Trays

- Tray holds 8 slides
- Pack of 4 trays and 1 lid

DESCRIPTION	CAT. NO.	SPECIAL
Disposable Slide Trays	HW-HS15951A	\$32.50



Microscope Slide Boxes®

- Cork interior
- Colour-coded for easy identification

DESCRIPTION	CAT. NO.	SPECIAL
25-Place Slide Box, Blue	HW-HS15989A	\$4.50
50-Place Slide Box, Blue	HW-HS15996A	5.00
100-Place Slide Box, Blue	HW-HS15994A	5.50

